Talk 4 and 8 Icon strips	Cut carefully along dotted lines, this sheet makes 3 icon strips
1. Before printing this template, make sure your page scaling is set to NONE and print orientation is set on LANDSCAPE. 2. Print or copy your loop sheet onto a piece of 8 1/2" by 11" sticky backed label pager.	

Print or copy your lcon sheet onto a piece of 8 1/2" by 11" sticky backed label paper.
 Cut out the lcons you need, peel off the backing and press it in place within the lcon template guidelines. Make sure the lcon square sticker edges are completely flat.
 Cut your lcon strip out along the dotted lines. When inserting your icon strip into the communicator hold down the leading edge so it will not get caught on the cell dividers. If needed, a thin laminate sheet such as those used to protect book covers may be applied over your icon strip to protect it.